easy fill

Bond LC

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Bond LC

Instructions for use

Light curing one component bonding agent for the total-etch-technique

PART 1: Introduction

Bond LC is a light curing, one part primer and bonding system which does not require additional mixing. Bonding is carried out in a single stage. Bond LC was developed for firmly bonding composites and compomers to enamel and dentine as well as non-precious and precious metals. Bond LC is also indicated for impregnating root canals prior to filling or cementing passive or threaded root posts.

Bond LC can also be used as a primer for indirect adhesive restorations, e.g. porcelain and composite inlays, onlays, laminate veneers, crowns and bridges placed with chemically or dual curing cement.

The high bond strength between Bond LC and enamel or dentine is based on the same principles as apply to glasionomer cements. The methacrylic carboxylic acid ester provides for long term high bond strength and biocompatibility. In case of irritation or methacrylate allergy, do not use this material.

Bond LC is compatible with all commonly available brands of light curing composite. As CercoBond is based on ethanol and hydrophilic, it can be applied to slightly moist dentine surfaces using the "Wet Bonding" technique.

In cases where a dual-curing or self-curing priming and bonding

system is recommended, Bond LC can be dual-cured using a 1 : 1 mixture of Bond LC and Bond LC Activator. Bond LC Activator is available separately. This special type of application is described in the Bond LC Activator instructions.

PART 2: Using Bond LC for filling materials

Preparing for using light curing composite and compomer filling materials

1.1 Isolating

It is advisable to use a rubber dam.

1.2 Cavity preparation

Before beginning preparation, clean the tooth with pumice and water. To ensure that as little healthy tooth structure as possible is reduced, only prepare the tooth minimally. Enamel cavity margins must be bevelled slightly (0.5 – 1.0 mm) to enlarge the adhesive surface in the enamel and increase the bond strength.

1.3 Pulp protection

In very deep cavities those areas in close proximity to the pulp must be coated with a thin layer of calcium hydroxide cavity liner.

1.4 Etching the enamel and dentine

Apply Etching Gel onto the enamel and dentine surfaces beginning with the enamel bevels. Leave the etching gel in place for 20 seconds. On primary teeth and teeth

high in fluoride, a 60 second etch is recommended. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream but do not desiccate. A slightly wet dentine surface is important for the function of Bond LC. The etched enamel bevel should have a chalky white appearance.

Precautions when etching

The etched enamel surface must not be touched or contaminated with saliva prior to applying the Bond LC. Should it become contaminated, repeat the procedure described above, i.e. etch, rinse with water and dry again.

Warning: Do not allow etching gel to contact the mucosa, eyes or skin. Should this occur inadvertently, rinse immediately with copious amounts of water.

Although it is usually not necessary to etch the enamel for typical compomer applications, etching the enamel does increase the bond strength considerably.

1.5 Applying Bond LC

Shake the Bond LC before use.

Use a brush to apply copious amounts of Bond LC to the dentine and enamel and brush it in vigorously for 30 seconds. Then dry

carefully with oil-free compressed air for approximately 15 seconds.

Cure the Bond LC coat for 20 seconds with a Halogen lamp before applying a second coat.

CAUTION: Do not spill Bond LC! If Bond LC is not used immediately, keep it in subdued light to prevent it curing prematurely due to ambient light. Bond LC will not cure spontaneously.

1.6 Applying the second coat of Bond LC

Use a brush to apply a second coat of Bond LC to the surfaces being adhered and brush it in vigorously for approximately

30 seconds. Then dry for approximately 15 seconds with oil-free

compressed air and cure for 20 seconds with a Halogen lamp.

The light curing composite filling can then be placed.

1.7 Placing, curing and finishing the filling

The CapoHybrid instructions describe how the composite filling is placed, cured and finished.

2. Priming the root canals

2.1 Preparation

It is essential to isolate the teeth with a rubber dam. Clean the teeth with hydrogen peroxide solution (3%) or sodium hypochlorite solution. When placing passive or threaded endodontic posts with

luting composites, refer to the manufacturer's instructions. Prepare the canal to determine the length of root post required. Adapt the post to the canal and reduce the length if necessary. Remove the post, keep it in a safe place and prepare the canal for cementing the post.

2.2 Priming the root canals

As the surfaces of the root canal cannot be reached with a brush, Bond LC must be injected into the prepared canals with a syringe to rinse them for 60 seconds in four separate stages. Ensure that all surfaces being bonded

are wetted continually. Wherever possible, use a brush to work the Bond LC in and enhance wetting. Soak up the excess Bond LC from the root canal using, for example, a soft paper point and dry the treated surfaces carefully for 15 seconds with oil-free compressed air. The Bond LC is then cured for 20 seconds with a Halogen lamp.

2.3 Applying the luting composite

Place and cure the luting composite as described in its manufacturer's instructions.

PART 3: Special applications for Bond LC

3. Using Bond LC on precious/non-precious metals and porcelain for repairing crown and bridgework

3.1 Conditioning the surfaces of the restoration

Roughen the surfaces of the metal or fractured porcelain with

a diamond instrument (40 mm) – if necessary, sandblast them in

addition (with an intraoral sandblaster).

Apply Etching to the rough surfaces, rinse them with copious

amounts of water and dry them with oil-free, dry compressed air.

Warning: Do not allow etching gel to contact the mucosa, eyes

or skin. Should this occur inadvertently, rinse immediately with

copious amounts of water.

3.2 Applying Bond LC

Brush Bond LC onto the metal surface carefully and dry it for

15 seconds with oil-free, dry compressed air. Cure the coat with a

dental curing lamp for 20 seconds. Begin repairing the restoration.

PART 4: Storage and technical data

4. Storage

Store at 10-25°C.

5. Technical data

Light curing time (with a curing lamp) 20 seconds

Chemical curing time Approx. 3 minutes

(1 : 1 mixture of Bond LC and Bond LC Activator, excluding light and air)

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und Innovationen GmbH